



**ACTUAL SIZE\***

\*On A4 Print

# QT DEVICE

The QT Device delivers a richer and more fulfilling response experience than ever before. Providing the same features as other RF clickers, the QT expands upon functionality by offering advanced capabilities such as a full QWERTY keypad and backlit screen for easy viewing.

## SPECS

### ENCLOSURE

- Dimensions (LxWxH): 143mm x 73mm x 19mm
- Weight: 136 grams (with batteries)
- Lightweight, compact format
- LCD screen with backlight for optimum viewing
- Durable, rugged case ensures maximum longevity under extreme conditions

### USER INPUT

- QWERTY Entry
- 56 Keys, Glow in the Dark Keypad
- Supported Symbols: . , ` " ; : ! \ ? @ # \$ % ^ & \* ( ) \_ + - = / < > [ ] { } ? € £

### DISPLAY

- Easy to read LCD Display with 57.1mm x 26.9mm viewable area
- Channel Number, Battery Life Indicator, Answer Choice, Response Confirmation

### POWER & POWER MANAGEMENT

- Two AAA batteries
- Average battery life is 6 to 12 months

### USER IDENTIFICATION

- Fixed device ID
- Use ResponseCard Programmer to adjust ID with Insight 360 Setup is plug-and-play with TurningPoint

### RF TECHNOLOGY

- Up to 74 sessions can be running at one time in close proximity without interference
- Fully FCC, CE and Industry Canada certified

### RANGE

- Presentation Mode - Range for one RF Receiver is 60m (120m coverage)
- Test Mode - Data is stored in non-volatile memory until individual finishes test and uploads responses

### SUPPORTED PRODUCTS

- TurningPoint v.5.3
- TurningPoint Cloud
- Insight 360 v.2.3.2 and higher
- Flow v.1.7 and higher

## FUNCTIONALITY HIGHLIGHTS

### INPUT CAPABILITIES

- Multiple Choice, Multiple Response, Numeric Response, Short Answer, True/False, Essay

### PARTICIPANT LIST CREATION

- Automatically create participant lists with QT Device's user stored data or have participants log in to an existing list.
- Easily move from one session to another.

### SELF-PACED TESTING

- Participants use QT Device to navigate to any question on the test or assignment and select a response. The QT Device has the largest storage capacity in the industry, storing all responses into non-volatile memory to ensure data security and stability.

### SUCCESSFUL TRANSMISSION

- Participants and presenters receive visual successful transmission notification.